

| | | | | |
|------------------------------------|-----------------------------------|--|---------------------------------------|---------------------------------|
| Алматы (7273)495-231 | Иваново (4932)77-34-06 | Магнитогорск (3519)55-03-13 | Пермь (342)205-81-47 | Тверь (4822)63-31-35 |
| Ангарск (3955)60-70-56 | Ижевск (3412)26-03-58 | Москва (495)268-04-70 | Ростов-на-Дону (863)308-18-15 | Тольятти (8482)63-91-07 |
| Архангельск (8182)63-90-72 | Иркутск (395)279-98-46 | Мурманск (8152)59-64-93 | Рязань (4912)46-61-64 | Томск (3822)98-41-53 |
| Астрахань (8512)99-46-04 | Казань (843)206-01-48 | Набережные Челны (8552)20-53-41 | Самара (846)206-03-16 | Тула (4872)33-79-87 |
| Барнаул (3852)73-04-60 | Калининград (4012)72-03-81 | Нижний Новгород (831)429-08-12 | Саранск (8342)22-96-24 | Тюмень (3452)66-21-18 |
| Белгород (4722)40-23-64 | Калуга (4842)92-23-67 | Новокузнецк (3843)20-46-81 | Санкт-Петербург (812)309-46-40 | Ульяновск (8422)24-23-59 |
| Благовещенск (4162)22-76-07 | Кемерово (3842)65-04-62 | Ноябрьск (3496)41-32-12 | Саратов (845)249-38-78 | Улан-Удэ (3012)59-97-51 |
| Брянск (4832)59-03-52 | Киров (8332)68-02-04 | Новосибирск (383)227-86-73 | Севастополь (8692)22-31-93 | Уфа (347)229-48-12 |
| Владивосток (423)249-28-31 | Коломна (4966)23-41-49 | Омск (3812)21-46-40 | Симферополь (3652)67-13-56 | Хабаровск (4212)92-98-04 |
| Владикавказ (8672)28-90-48 | Кострома (4942)77-07-48 | Орел (4862)44-53-42 | Смоленск (4812)29-41-54 | Чебоксары (8352)28-53-07 |
| Владимир (4922)49-43-18 | Краснодар (861)203-40-90 | Оренбург (3532)37-68-04 | Сочи (862)225-72-31 | Челябинск (351)202-03-61 |
| Волгоград (844)278-03-48 | Красноярск (391)204-63-61 | Пенза (8412)22-31-16 | Ставрополь (8652)20-65-13 | Череповец (8202)49-02-64 |
| Вологда (8172)26-41-59 | Курск (4712)77-13-04 | Петрозаводск (8142)55-98-37 | Сургут (3462)77-98-35 | Чита (3022)38-34-83 |
| Воронеж (473)204-51-73 | Курган (3522)50-90-47 | Псков (8112)59-10-37 | Сыктывкар (8212)25-95-17 | Якутск (4112)23-90-97 |
| Екатеринбург (343)384-55-89 | Липецк (4742)52-20-81 | | Тамбов (4752)50-40-97 | Ярославль (4852)69-52-93 |

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.bbBraun.nt-rt.ru | | bng@nt-rt.ru

Технические характеристики на прочие инструменты и принадлежности компании **B. BRAUN**

Виды товаров: наборы катетеров для диагностического перитонеального лаважа, желудочные зонды, назогастральные и назоинтестинальные трубки из ПВХ, ректальные катетеры, ирригационные катетеры, закрытые системы для повторного лаважа и диагностики, наборы для оршения для послеоперационной ирригации, наборы удлинителей, наборы для инсуффляции кислорода, наборы для аспирации, вакуумные аспирационные камеры для сбора жидкостей, отсасывающие катетеры, устройства для отсасывания слизи, устройства для отсасывания слизистой для новорожденных и др.

Catheter set for diagnostic peritoneal lavage



Peritofix

Set for diagnostic peritoneal lavage consisting of:

Cutfix® scalpel fig. E 11

Puncture trocar Ch. 10

- With adjustable insertion length
- Metal trocar
- Plastic cannula

Certon® catheter, Ch. 9

- Made of polyurethane, length 50 cm
- Closed tip, lateral openings
- Coloured length markings

Gastric Tubes

Naso-gastric and naso-intestinal PVC tubes for short term use only.



Ernaehrungssonde_4396405_



Feeding Tubes

For **short term** application of tube feeding diets into the stomach.

Features

- Radiopaque (stripe)
- Luer (L) & tunnel (TR) connector
- Gastric tube CH16 equipped with stylet (Mandrin M)

Rectal catheter

Irrigation catheters



Rectal Catheter

- Color-coded cone fitting
- Closed rounded tip, 2 lateral eyes

Gastric Tube with stylet

PVC tube



Stomach Tube

Gastric PVC tube for short term use.

Features

- Radiopaque (stripe)
- With stylet
- Funnel connector

Rectobag®

Closed system for repeated lavage and diagnostics



Rectobag_



Rectobag

Closed system for repeated intestinal lavage and diagnostics

- Conical connector
- Large-bore transfer tube (ø 7.5 mm), length 150 cm with closing clamp
- Lavage bag 2 l with graduation
- Liquid and odor-tight inlet
- Sturdy suspension eyelets

Glycine 1.5% B. Braun

Irrigation solution



Glycine 1,5% B. Braun

Sterile and pyrogen-free

Advantages

Ecobag® - the container with top performance

- Closed system for unvented use
- Needle-free product system
- 2 sterile, resealing ports
- PVC/latex-free
- Unbreakable

Features

Presentations:

- Ecobag® CE PVC-free bag 3000 ml, 5000 ml
- Ecobag® click PVC-free bag 3000 ml, 5000 ml
(may differ from country to country)

Indication

- Intra- and postoperative applications include irrigations which are performed during endoscopic examinations and operations in body cavities and joints, conducted in the fields of Urology, Orthopedics, Accident Surgery, and Gynecology.

Ringer B. Braun

Irrigation solution



Ringer B. Braun



Ringer B. Braun Irrigation Solution



Sterile and pyrogen-free

Advantages

Ecobag® - the container with top performance

- Closed system for unvented use
- Needle-free product system
- 2 sterile, resealing ports
- PVC/latex-free
- Unbreakable

Features

Presentations:

- Ecobag® CE PVC-free bag 2000 ml, 3000 ml and 5000 ml
- Ecobag® click PVC-free bag 2000 ml, 3000 ml and 5000 ml
(may differ from country to country)

Indication

Intra- and postoperative applications include irrigations which are performed during endoscopic examinations and operations in body cavities and joints, conducted in the fields of Urology, Orthopedics, Accident Surgery, and Gynecology.

Ecoclick® No. 1

Irrigation Set



Ecoclick No. 1

Single line irrigation set for intra- and postoperative irrigation of the bladder with Ecobag® click.

Advantages

Ecoclick® can be connected spike-free and with a simple click to the Ecobag®click green connector. Once connected, the Ecoclick® tubing rotates on its own axes for flexibility during handling.

Features

- Ecoclick®-connector
- Roller clamp
- Drip chamber

- Standard cone with 15 cm silicone tubing for connection with resectoscopes
- Flow rate approx. 1000 ml/min
- Length approx. 210 cm

Indication

Intra- and postoperative urological interventions.

Ecoclick® No. 4

Irrigation set



Ecoclick No. 4

Single line irrigation set with non return valve for intra- and postoperative urological irrigation with Ecobag® click

Advantages

Ecoclick® can be connected spike-free and with a simple click to the Ecobag®click green connector. Once connected, the Ecoclick® tubing rotates on its own axes for flexibility during handling.

Features

- Ecoclick® connector
- Drip chamber
- Roller clamp

- Non-return valve
- Standard cone with 15 cm silicone tubing for connection with resectoscopes
- Flow rate approx. 1000 ml/min
- Length approx. 210 cm

Indication

Intra- and postoperative urological interventions.

Ecoclick® Nr. 2 Arthroscopy

Irrigation set



Ecoclick - Arthros No. 2

Single line irrigation set for arthroscopical use with Ecobag® click

Advantages

Ecoclick® can be connected spike-free and with a simple click to the Ecobag®click green connector. Once connected, the Ecoclick® tubing rotates on its own axes for flexibility during handling.

Features

- Ecoclick® connector
- Drip chamber
- Roller clamp
- Luer-Lock fitting for connection with arthroscopes
- Length approx. 330 cm

Indication

Arthroscopy

Ecospike Nr. 1

Irrigation set



Ecospike No. 1

Single line irrigation set for intra- and postoperative urological irrigation

Features

- Drip chamber
- Roller clamp
- Standard cone with 15 cm silicone tubing for connection with resectoscopes
- Flow rate approx. 800 ml/min
- Tube length approx. 210 cm

Indication

Intra- and postoperative urological interventions.

Ecospike Nr. 4

Irrigation Set



Ecospike No. 4

Single line irrigation set with non return valve for intra- and postoperative urological irrigation.

Features

- Drip chamber
- Roller clamp
- Non-return valve
- Standard cone with 15 cm Silicone tubing for connection with resectoscopes
- Flow rate: approx. 800 ml/min
- Tubing length till Y-connection: 50 cm each
- Length approx. 210 cm

Indication

Intra- and postoperative urological interventions.

NaCl 0.9% B. Braun

Irrigation solution



NaCL 0.9 B. Braun

Sterile and pyrogen-free

Advantages

Ecobag® - the container with top performance

- Closed system for unvented use
- Needle-free product system
- 2 sterile, resealing ports
- PVC/latex-free
- Unbreakable

Features

Presentations:

- Ecobag[®] CE PVC-free bag 2000 ml, 3000 ml and 5000 ml
- Ecobag[®] click PVC-free bag 2000 ml, 3000 ml and 5000 ml
(may differ from country to country)

Indication

Intra- and postoperative applications include irrigations which are performed during endoscopic examinations and operations in body cavities and joints, conducted in the fields of Urology, Orthopedics, Accident Surgery, and Gynecology.

Aqua B. Braun

Intraoperative irrigation solution



Aqua B. Braun in Ecobag

Advantages

Ecobag® - the container with top performance

- Closed system for unvented use
- Needle-free product system
- 2 sterile, resealing ports
- PVC/latex-free
- Unbreakable

Features

Presentations:

- Ecobag® 2000 ml and 3000 ml
(may differ from country to country)

Indication

Intra- and postoperative applications include irrigations which are performed during endoscopic examinations and operations in body cavities and joints, conducted in the fields of Urology, Orthopedics, Accident Surgery, and Gynecology.

Ecoclick® No. 2

Irrigation set



Ecoclick No. 2

Twin line irrigation set for intra- and postoperative irrigation of the bladder with Ecobag® click.

Advantages

Ecoclick® can be connected spike-free and with a simple click to the Ecobag®click green connector. Once connected, the Ecoclick® tubing rotates on its own axes for flexibility during handling.

Features

- Ecoclick® -connectors
- 2 pressure clamps
- Y-connection

- Drip chamber
- Roller clamp
- Standard cone with 15 cm silicone tubing for connection with resectoscopes
- Flow rate approx. 1000 ml/min
- Tubing length till Y-connection approx. 50 cm each
- Length approx. 210 cm

Indication

Intra- and postoperative urological interventions.

Ecoclick® Nr.1 Arthroscopy

Irrigation set



Ecoclick - Arthros No. 1

Single line irrigation set for arthroscopical use with Ecobag® click

Advantages

Ecoclick® can be connected spike-free and with a simple click to the Ecobag®click green connector. Once connected, the Ecoclick® tubing rotates on its own axes for flexibility during handling.

Features

- Ecoclick® connector
- Drip chamber
- Roller clamp
- 15 cm standard silicone tubing
- Length approx. 330 cm

Indication



Arthroscopy

Extension set for Ecospike/Ecoclick®

Extension set with Luer-Lock connector



Fitting set Luer Lock



Fitting set Luer-Lock

The adapter fits to standard silicone tubing sets of Ecospike and Ecoclick® sets and provides a Luer-Lock connection to the surgical instrument.

Features

- Cone fitting
- Luer-Lock fitting
- Length approx. 60 cm

Ecospike Nr. 2

Irrigation Set



Ecospike No. 2

Twin line irrigation set for intra- and postoperative urological irrigation.

Features

- 2 Ecospike connectors
- 2 pressure clamps
- Y-connection
- Drip chamber
- Roller clamp

- Standard cone with 15 cm Silicone tubing for connection with resectoscopes
- Flow rate: approx. 800 ml/min
- Tubing length till Y-connection: 50 cm each
- Length approx. 210 cm

Indication

Intra- and postoperative urological interventions

Ecospike Nr. 2 Arthroscopy

Irrigation Set



Ecospike-Arthros No. 2

Single line irrigation set for arthroscopical use.

Features

- Drip chamber
- Roller clamp
- Luer-Lock fitting for connection with arthroscopes
- Length approx. 330 cm

Indication

Ringer Lactate

Irrigation solution



Ringer Lactate

Sterile and pyrogen-free

Advantages

Ecolav® - easy to handle and time saving

- The soft container material and the spray nozzle allow precise pressure irrigation immediately after opening the "twist-off" cap
- No more troublesome, time consuming decanting of irrigation solutions and aspirating the solution into a syringe
- Precise irrigation with the spray nozzle

- Immediate pressure irrigation due to soft container material

Features

Presentations:

- Ecolav[®] PE bottle 30 ml, 100 ml and 500 ml
(may differ from country to country)

Indication

- Moistening of tissue and cleansing of the operating area during surgical procedures
- Rinsing of wounds in the event of injuries and burns
- Cleansing during stoma care
- Moistening of wound tamponades, cloths, bandages, and dressings

Ringer B. Braun

Irrigation Therapy



Ringer B. Braun

- Irrigation Solution
- Sterile and pyrogen-free

Aqua B. Braun

Topical irrigation solution



Aqua B. Braun

Sterile and pyrogen-free

Advantages

Ecolav® - easy to handle and time saving

- The soft container material and the spray nozzle allow precise pressure irrigation
- Immediately after opening the "twist-off" cap
- No more troublesome, time consuming decanting of irrigation solutions and aspirating the solution into a syringe
- Precise irrigation with the spray nozzle

- Immediate pressure irrigation due to soft container material

Features

Presentations:

- Ecolav[®] PE bottle 100 ml, 250 ml, 500 ml and 1000 ml
(may differ from country to country)

Indication

- Moistening of tissue and cleansing of the operating area during surgical procedures
- Rinsing of wounds in the event of injuries and burns
- Cleansing during stoma care
- Moistening of wound tamponades, cloths, bandages, and dressings
- Cleansing of operating instruments and accessories
- Filling of humidifiers and vaporizers

Oxygen Insufflation Set



Oxygen Insufflation Set

- Made of soft DEHP-free PVC
- Atraumatic nasal tip

Advantages

- Soft and short nasal connector minimally irritates mucosa and is comfortable for the patient
- Soft funnel enables easy connection to oxygen source

Vacufix® Suction Set

Vacuum suction system for collection of liquids



Vacufix® Suction Sets

Vacufix® Suction set consisting of:
Handle and tubing

- Profile and integrated finger grip for ergonomic touch
- Comfortable finger support
- Appropriate tip for each application
- Various connectors
- Safe connection
- Compatible with all common suction systems

Vacufix® Chamber

Vacuum suction system for collection of fluids



Vacufix® Chamber

- Provides a constant vacuum level
- Housing made of impact-resistant polyurethane
- Glass pane made of Securit
- Practical handle for secure handling
- Pressure stabilised to - 0.9 bar
- Disinfectant-resistant

Vacufix® Gelling Granulate 2 kg



Vacufix® Gelling Granulate 2 kg

Absorbs suctioned fluids prior to disposal of the Vacufix® collection bag

- With measuring cup

Suction Catheter Invasoft



Invasoft Suction Catheter

Straight suction catheter with air-vent tip, central opening, 2 lateral eyes

- Soft tubing of medical quality
- Easy gliding characteristics
- Smooth tubing surface and atraumatic rounded catheter tip for minimal mucous membrane irritation
- Air-vent tip combined with small lateral eyes prevents adhesion of the catheter tip
- Colour-coded cone fitting
- Soft PVC 75 ± 3 shore A

Indication

Suction catheters are flexible, long tubes used to remove respiratory secretions from the airway.

Suction Catheter Ideal

Flexible, long tube to remove respiratory secretions from the airway



Ideal Suction Catheter

Soft tubing of medical quality

- Easy gliding characteristics
- Smooth tubing surface and atraumatic rounded catheter tip for minimal mucous membrane irritation
- Small lateral eyes prevent adhesion of the catheter tip
- Colour-coded cone fitting
- Soft PVC 71 \pm 3 shore A

Suction Catheter Standard



Suction Cathetr Standard

- Tubing of medical quality
- Easy gliding characteristics
- Smooth tubing surface and atraumatic rounded catheter tip for minimal mucous membrane irritation
- Color-coded cone fitting
- Soft PVC (80 ± 3 Shore A)

Indication

Suction catheters are flexible, long tubes used to remove respiratory secretions from the airway.

Mucex

Mucous suction device for newborns



Mucex

For neonatology:

- Volume 25 ml
- Graduated

Finger tip



FingerTip

Finger tip for intermittent suction

Trachea Secretion Sampling Set

Trachea-secretion test set is used for obtaining tracheal secretion for analysis

- Breakproof sampling tube, 10 ml volume
- Patient data label
- Integrated finger tip
- Cone fitting Ch. 14
- Separate closing stopper

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.bb Braun.nt-rt.ru | | bng@nt-rt.ru